

## APPENDIX B

### “Finding the Characteristic Times, $T_1$ and $T_2$ ”

Source: Rasa Rejali (2014)



# Finding the Characteristic Times, $T_1$ and $T_2$

July 2014

## Determining $T_1$

In order to measure the spin lattice relaxation time, a  $180^\circ$  A pulse must be applied followed by a  $90^\circ$  B pulse, a time  $\tau$  later. You will only use two pulses, rather than a sequence of them. Note that the spectrometer must be properly tuned and the temperature of the magnet must be perfectly stable to achieve reliable results.

1. Ensure the spectrometer has been tuned according to the manual instructions; take note of the final values for the  $90^\circ$  pulse length and  $180^\circ$  pulse length found during the calibration process.
2. Set A-len to the  $180^\circ$  pulse time; set B-len to the  $90^\circ$  pulse time
3. You will vary  $\tau$  throughout the measurement of  $T_1$ , but a good starting point can be found in Table 1 for a few selected samples. However, keep in mind that these starting points are simply guidelines.
4. Adjust the period (P) according to the values in Table 1, when relevant. Otherwise, try modifying this parameter until you observe the appropriate signal.

Table 1: Data collection:  $\tau$  range and period for selected samples

Sample	Start value for $T_1$	Start value for $T_2$	Period
Light mineral oil	40 ms	10 ms	1.0 s
Heavy mineral oil	20 ms	10 ms	1.0 s
Pure water	1.4s	10 ms	5.0s (for $T_1$ ), 1.0s (for $T_2$ )

5. Set the number of B-pulses in a sequence (N) to 1.
6. Set the Sync toggle switch to B, and turn on the B Pulse toggle switch. You should observe a signal at this point; adjust the oscilloscope window as needed.




Figure 1: Observing the FID amplitude of the  $180^\circ$  pulse, following the  $90^\circ$  pulse a time  $\tau = 0.1s$  later.

- Using the Measure function on the oscilloscope, measure the peak after the  $90^\circ$  pulse.
- Continue measuring the peak of the  $90^\circ$  pulse as you vary  $\tau$ . Plot the FID amplitude as a function of  $\tau$  to determine  $T_1$ . For further information, please refer to Appendix A, pg 7-8.

### Determining $T_2$

To measure the spin-spin relaxation time, a  $90^\circ$  A pulse must be applied followed by a  $180^\circ$  B pulse, a time  $\tau$  later; an echo is observed a time  $\tau$  after the  $180^\circ$  pulse. By applying a sequence of B pulses, the characteristic time  $T_2$  can be determined. It may be useful to think of  $T_2$  as the time it takes for the FID echo amplitude to drop to 37%.

- Set A-len to the  $90^\circ$  pulse; set B-len to the  $180^\circ$  pulse; adjust the period (P) according to Table 1.
- Start with one B pulse in the sequence, and ensure the output pulse corresponds to Figure 2. You can then begin to increase the number of pulses in the sequence until you observe an exponential decay of the echo amplitude (refer to Figure 3).




Figure 2: A  $90^\circ$  pulse, followed by a  $180^\circ$  pulse, followed by an echo; each peak is separated by a time  $\tau = 0.1s$ .




Figure 3: Carr-Purcell principle of multiple echoes, exhibiting exponential decay of the echo amplitude.

- Ensure to apply the Meiboom-Gill correction to the Carr-Purcell sequence, such that the uncertainty does not propagate so as to create an attenuated echo.
- Please refer to Appendix A, pg 4-7 for further instruction.